



Sample Submission and Results Reporting

CONTACT INFORMATION

New Customer

Returning Customer

Sample Contact Person: _____

Account Name: _____ Account Number: _____

Phone: _____

Results will be emailed to the address provided and will be stored on your maxill online account. Login to access your results at any time.

Email sample results to (Name): _____

Email: _____

Call your H₂O Customer Service Representative 24 HOURS BEFORE taking water samples.
Dawn Giesbrecht 1-800-268-8633 ext. 240 | dawng@maxill.com

Instructions (Please read ALL instructions BEFORE taking water samples)

To help ensure accurate test results, your maxill H₂O Customer Service Representative will help co-ordinate your water sampling and submission. For accurate results, water samples must be kept cold at all times and received by maxill within 24 hours of sampling. Water samples can only be taken Monday or Tuesday and must be received by the maxill lab by 3 PM on Wednesday.

Pre-sampling Instructions


- Once you receive your Dental Unit Waterline (DUWL) Test Kit remove the cooler from the shipping box. Do not discard the shipping box, it will be used to return the cooler and samples for testing.
- Remove the ice pack from the cooler and fill it approximately 3/4 full with water. Lay the ice pack as flat as possible into the bottom of the cooler and place the cooler into the freezer, sitting flat. The ice pack should be allowed at least 24 hours to freeze completely.
- Before taking your water samples please contact Dawn**, your maxill H₂O Customer Service Specialist at 1-800-268-8633, extension 240 to co-ordinate when you will be taking your samples so your Priority Return pick-up time can be scheduled.
- Fill out the sample identification fields on the reverse of this page and on the test tube labels provided with your kit. Test tube labels are to be applied at time of sampling. Reference your maxill H₂O DUWL Cleaning Log for the correct area and source codes to ensure accurate logging. A digital copy of the cleaning log can be downloaded on the H₂O tab of your maxill online account.

Sampling Instructions

- Wearing examination gloves that have been disinfected (we recommend 70% isopropyl alcohol), remove hand-pieces from dental water line tubing.
- Disinfect the end of the tubing with an appropriate disinfectant (we recommend 70% isopropyl alcohol).
- Allow water to be flushed through the tubing for 20-30 seconds. Being careful to not touch or breath on the interior of the sample tube, fill the tube to approximately 3 cm below the top of the tube.
- Firmly apply the cap and apply the appropriate label to the lower half of the test tube with the information reading lengthwise on the tube (fig. 1).
- Place tubes into the provided foam envelope and refrigerate at 1°-8°C until ready to ship (Ideally, samples should be taken within 1-2 hours of shipping).

Packing and Shipping instructions

- Before shipping, gently remove the frozen ice pack from the cooler.
- Place foam envelope containing the test tubes into the cooler. Place the ice pack on top of the test tubes (fig. 2).
- Attach the lid on the cooler and place the cooler into the provided shipping box.
- Seal the box and attach the provided Priority Return shipping label.



PROPOSED SAMPLE DATE:

DATE SAMPLED:

| |
|----------------|
| YYYY / MM / DD |
|----------------|

| Sample Identification | | | FOR maxill USE ONLY - Total Viable Count Testing | | |
|-----------------------|----------------|-------------|--|--------|-------------|
| # | OP / Area Code | Source Code | ID | Result | Action Code |
| 1 | | | | | |
| 2 | | | | | |
| 3 | | | | | |
| 4 | | | | | |
| 5 | | | | | |
| 6 | | | | | |
| 7 | | | | | |
| 8 | | | | | |
| 9 | | | | | |
| 10 | | | | | |
| 11 | | | | | |
| 12 | | | | | |
| 13 | | | | | |
| 14 | | | | | |
| 15 | | | | | |
| 16 | | | | | |
| 17 | | | | | |
| 18 | | | | | |
| 19 | | | | | |
| 20 | | | | | |
| 21 | | | | | |
| 22 | | | | | |
| 23 | | | | | |
| 24 | | | | | |
| 25 | | | | | |
| 26 | | | | | |
| 27 | | | | | |
| 28 | | | | | |
| 29 | | | | | |
| 30 | | | | | |

*Failure to include all information may result in longer wait times. It is imperative that samples be kept cold in transit to maxill.
TNTC = Too numerous to count.
In the event a result cannot be obtained due to the number of colonies.

| | |
|----------------------------|--|
| Date Received: | Tested By: |
| Test Date: | Flag Results: <input type="radio"/> Yes <input type="radio"/> No |
| Results Entered By: | |